


Fitting a cylinder to path control points using the SimVascular Python API

This document describes how to use the SimVascular Python API to fit a cylinder to two control points defined on a path using the **fit_cylinder_to_path.py** Python script.

The **fit_cylinder_to_path.py** script defines a Python class named **PathFitCyl** that encapsulates the data and functions used to create a cylinder aligned with the axis defined by two control points on a path.

A path with control points marked by blue squares. Control points are indexed using integer IDs: 0,1,2,...,26.

Creating a cylinder between controls points 7 (red) and 8 will have an axis aligned with those points (red arrow) and length equal to the distance between those points.



Executing the **fit_cylinder_to_path.py** script defines the **PathFitCyl** class. This class can then be used to create **PathFitCyl** objects. An **PathFitCyl** object is created by given a path name to the object constructor

`object = PathFitCyl(path_name)`, where *path_name* is the name of an SV path.


A cylinder is created using the **fit()** method supplying two control points integer IDs and the cylinder radius

`object.fit(id1, id2, radius)`

This creates a cylinder and displays it in the SV graphics window.


The following example uses the SimVascular Demo project.

Open Python console



We will use the **aorta** path (red arrow). There are 27 controls points indexed 0,1,2,...26. Control points are shown as blue squares.

Open the Python console by selecting the Python icon located on the tool bar in the upper right of the SV window.




Select Text Editor


We now want to read in the the **fit_cylinder_to_path.py** Python script.

Select the **Text Editor** button at the bottom of the SV Python Console.

Select to read in a Python script




Select the  icon to open a file browser used to select the script.



Select Python the script to read in.


Select the play button to execute the script.



Select the icon to execute the script which defines the **PathFitCyl** class.

Select **Console** to go back to the Python console.

Select the **Console** button to go back to the Python console.



To create a cylinder type the follow commands into the Python console


```
>>> aorta_fit = PathFitCyl('aorta')
[PathFitCyl] Repository name: aortaRepo
[PathFitCyl] Writing to default directory /
[PathFitCyl] **** ERROR: Can't write to directory: /
[PathFitCyl]           Set directory using the 'set_write_directory' method.
[PathFitCyl] Add aortaRepo to the repository.
[PathFitCyl] Number of control points: 27
[PathFitCyl] Number of path points: 616
>>> aorta_fit.set_write_directory('/Users/parkerda/')
>>> aorta_fit.fit(7,8,1.0)
[PathFitCyl] Distance between control points 7 8 : 2.38058
[PathFitCyl] Point1: [-0.138486 1.017057 0.335150]
[PathFitCyl] Point2: [-0.241842 1.740550 -1.930475]
[PathFitCyl] Axis: [0.043416 -0.303914 0.951710]
[PathFitCyl] Center: [-0.190164 1.378804 -0.797662]
[PathFitCyl] Creating cylinder aorta_fit_cylinder_7_8
```

Create PathFitCyl object.


Set the directory used to write cylinder geometry.

The **fit_cylinder_to_path.py** script writes out the cylinder geometry to a VTK .vtp file and then reads it back in order to display it. This is a work around for a missing API function to create geometry directly.

The directory where the geometry files are written should be set after creating the **PathFitCyl** object using the **set_write_directory()** method. This directory must have write permissions for the current user.



This creates a cylinder named **aorta_fit_cylinder_7_8** under the **Repository** node of **SV Data Manager** (red arrow). The cylinder is displayed in the SV graphics window. Note that the cylinder is not immediately displayed, the graphics window but be zoomed or rotated for the geometry to appear.




Another cylinder can be created using

```

>>> aorta_fit.fit(15,16,1.0)
[PathFitCyl] Distance between control points 15 16 : 2.31655
[PathFitCyl] Point1: [-0.825675 0.186380 -10.954674]
[PathFitCyl] Point2: [-2.214444 -1.282020 -12.086725]
[PathFitCyl] Axis: [0.599498 0.633873 0.488679]
[PathFitCyl] Center: [-1.520059 -0.547820 -11.520700]
[PathFitCyl] Creating cylinder aorta_fit_cylinder_15_16

```




To modify the cylinder radius for control points (7,8) delete the **aorta_fit_cylinder_7_8**

Repository node of **SV Data Manager** (red arrow) by right clicking on **aorta_fit_cylinder_7_8** and selecting the **Remove** menu item.

You can now create another cylinder with a radius of 0.5

```
>>> aorta_fit.fit(7,8,0.5)
[PathFitCyl] Distance between control points 7 8 : 2.38058
[PathFitCyl] Point1: [-0.138486 1.017057 0.335150]
[PathFitCyl] Point2: [-0.241842 1.740550 -1.930475]
[PathFitCyl] Axis: [0.043416 -0.303914 0.951710]
[PathFitCyl] Center: [-0.190164 1.378804 -0.797662]
[PathFitCyl] Creating cylinder aorta_fit_cylinder_7_8
[PathFitCyl] Remove cylinder aorta_fit_cylinder_7_8
```



Position: <-0.79, 0.09, -10.99> mm; Index: <251, 26, 153>; Time: 0.00 ms; Pixel value: 163.00 5.20 GB (32.48 %)